

Anti-Rxrα Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Chip, Chip-Seq, Elisa - P19793

Item Number: CM0003365



Introduction

Recombinant rabbit monoclonal antibody against human RXRα (NR2B1), validated in WB, IP, IF/ICC, ChIP, ChIP-seq, and ELISA. Detects endogenous levels of RXRα (51 kDa) in human, mouse, and rat samples. Ideal for studying retinoid receptor signaling, gene regulation, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	RXRα Rabbit mAb
Remarks/Alias	NR2B1; RXRalpha; RXR-alpha; RXRα
Species	Human
GeneID (Human)	6256
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human RXRα (P19793)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, ChIP, ChIP-seq, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19793
Protein Weight	51 kDa
Shipping	Ice bag