

Anti-Rxrα Rabbit Polyclonal Antibody For Wb, Elisa - P19793

Item Number: CM0014999



Introduction

Rabbit polyclonal antibody against human RXRα (NR2B1). Validated for Western blot and ELISA, with reactivity to human, mouse, and rat. Ideal for transcriptional regulation and nuclear receptor studies.

[Learn More](#)

Parameter	Value
Product Name	RXRα Rabbit pAb
Remarks/Alias	NR2B1; RXRalpha; RXR-alpha; RXRα
Species	Human
Gene ID (Human)	6256
Gene ID	6256
Immunogen	Synthetic peptide within amino acids 1-100 of human RXRα (NP_002948.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19793
Protein Weight	51kDa
Shipping	Ice bag