


Anti-Rxrb Rabbit Monoclonal Antibody For Wb, Ip, Chip, Elisa - P28702

Item Number: CM0011342



Introduction

Rabbit monoclonal antibody targeting human RXRB (RXR-beta). Suitable for Western blot, immunoprecipitation, ChIP, and ELISA. Cross-reacts with mouse. Ideal for studying retinoic acid signaling and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	RXRB Rabbit mAb
Remarks/Alias	NR2B2, DAUDI6, RCoR-1, RXRbeta, H-2RIIBP, RXR-beta
Species	Human
Gene ID (Human)	6257
Gene ID	6257
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RXRB (NP_068811.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ChIP, ELISA
Cross Reactivity	Mouse
SWISS	P28702
Protein Weight	57 kDa
Shipping	Ice bag