


Anti-Rxrb Rabbit Polyclonal Antibody For Wb, Elisa - P28702

Item Number: CM0021587



Introduction

High-quality rabbit polyclonal antibody against human RXRB (RXR-beta, NR2B2), suitable for Western blotting and ELISA. Detects human, mouse, and rat RXRB. Ideal for studying retinoic acid signaling and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	RXRB Rabbit pAb
Remarks/Alias	NR2B2; DAUDI6; RCoR-1; RXRbeta; H-2RIIBP; RXR-beta; RXRB
Species	Human
GeneID (Human)	6257
GeneID	6257
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28702
Protein Weight	57kDa
Shipping	Ice bag