

## Anti-Ryr1 Polyclonal Antibody For Wb, Elisa - P21817

Item Number: CM0020837



### Introduction

Rabbit polyclonal antibody against human Ryanodine receptor 1 (RYR1). Validated for Western blot (WB) and ELISA. Cross-reacts with mouse. Suitable for muscle contraction, calcium signaling, and developmental biology research.

[Learn More](#)

| Parameter        | Value                                                                                                   |
|------------------|---------------------------------------------------------------------------------------------------------|
| Product Name     | RYR1 Rabbit pAb                                                                                         |
| Remarks/Alias    | CCO; KDS; MHS; RYR; MHS1; RYDR; SKRR; RYR-1; CMYP1A; CMYP1B; PPP1R137; RYR1                             |
| Species          | Human                                                                                                   |
| GeneID           | 6261                                                                                                    |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 800-900 of human RYR1 (NP_000531.2). |
| Source           | Rabbit                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                   |
| Application      | WB, ELISA                                                                                               |
| Cross Reactivity | Mouse                                                                                                   |
| SWISS            | P21817                                                                                                  |
| Protein Weight   | 565kDa                                                                                                  |
| Shipping         | Ice bag                                                                                                 |