


Anti-Ryr2 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q92736

Item Number: CM0005542



Introduction

Recombinant rabbit monoclonal antibody targeting ryanodine receptor 2 (RYR2), validated in WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for cardiac calcium signaling and arrhythmia research.

[Learn More](#)

Parameter	Value
Product Name	RYR2 Rabbit mAb
Remarks/Alias	RyR; ARVC2; ARVD2; RYR-2; VTSIP; VACRDS
Species	Human
Gene ID	6262
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92736
Protein Weight	565kDa
Shipping	Ice bag