


Anti-Ryr2 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q92736

Item Number: CM0016359



Introduction

High-quality rabbit polyclonal antibody targeting human RYR2 (UniProt Q92736). Validated for WB, IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for cardiac calcium signaling and arrhythmia research.

[Learn More](#)

Parameter	Value
Product Name	RYR2 Rabbit pAb
Remarks/Alias	RyR; ARVC2; ARVD2; RYR-2; VTSIP; VACRDS; RYR2
Species	Human
GeneID (Human)	6262
GeneID	6262
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 4868-4967 of human RYR2 (NP_001026.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
UniProt Accession	Q92736
Protein Weight	565 kDa
Shipping	Ice bag