

Anti-S100A1 Monoclonal Antibody For Wb, Ihc-P, Elisa - P23297

Item Number: CM0006908



Introduction

Rabbit anti-human S100A1 monoclonal antibody validated for WB, IHC-P, ELISA. Detects S100A1 in human, mouse, and rat samples.

Recombinant immunogen, 11 kDa calcium-binding protein ideal for cardiac and calcium signaling research.

[Learn More](#)

Parameter	Value
Product Name	S100/S100A1 Rabbit mAb
Remarks/Alias	S100; S100A; S100-alpha; S100/S100A1
Species	Human
Gene ID (Human)	6271
Immunogen	Recombinant protein of human S100/S100A1
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23297
Protein Weight	11kDa
Shipping	Ice bag