

Anti-S100A10 Monoclonal Antibody For Wb, Ihc-P, Elisa - P60903

Item Number: CM0003642



Introduction

Rabbit monoclonal antibody targeting human S100A10 (p11, calpactin I light chain). Validated for WB, IHC-P, and ELISA applications. An essential tool for investigating ANXA2 dimerization and phosphorylation regulation.

[Learn More](#)

Parameter	Value
Product Name	S100A10 Rabbit mAb
Remarks/Alias	42C; P11; p10; GP11; ANX2L; CAL1L; CLP11; Ca[1]; ANX2LG; S100A10
Species	Human
GeneID (human)	6281
GeneID	6281
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-97 of human S100A10 (P60903).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P60903
Protein Weight	11kDa
Shipping	Ice bag