

Anti-S100A11 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P31949

Item Number: CM0002486

Introduction

Rabbit recombinant monoclonal antibody to human S100A11, suitable for Western blot, immunohistochemistry, immunofluorescence, and ELISA applications. Designed for research on keratinocyte differentiation and cornification processes.

[Learn More](#)

Parameter	Value
Product Name	S100A11 Rabbit mAb
Remarks/Alias	MLN70; S100C; HEL-S-43; S100A11
Species	Human
GeneID (human)	6282
GeneID	6282
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P31949
Protein Weight	12 kDa
Shipping	Ice bag