


Anti-S100A11 Polyclonal Antibody For Wb, If/Icc, Elisa - P31949

Item Number: CM0021443



Introduction

Rabbit polyclonal antibody targeting human S100A11 (calgizzarin). Validated for WB, IF/ICC, ELISA. Suitable for keratinocyte differentiation and cornification studies. Reliable reagent for IVD and research applications.

[Learn More](#)

Parameter	Value
Product Name	S100A11 Rabbit pAb
Remarks/Alias	MLN70; S100C; HEL-S-43; S100A11
Species	Human
GeneID (human)	6282
GeneID	6282
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human S100A11 (NP_005611.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P31949
Protein Weight	12kDa
Shipping	Ice bag