


Anti-S100A11 Rabbit Polyclonal Antibody [KO Validated] For Wb, If/Icc, Ip, Elisa - P31949

Item Number: CM0025976



Introduction

High-quality anti-S100A11 rabbit polyclonal antibody validated for knockout (KO) experiments. Suitable for WB, IF/ICC, IP, ELISA in human samples. Targets calgizzarin involved in keratinocyte differentiation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] S100A11 Rabbit pAb
Remarks/Alias	MLN70; S100C; HEL-S-43; 11
Species	Human
GeneID	6282
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-105 of human S100A11 (NP_005611.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	P31949
Protein Weight	12kDa
Shipping	Ice bag