


Anti-S100A12 Polyclonal Antibody For IHC-P, Elisa - P80511

Item Number: CM0017730



Introduction

Rabbit polyclonal antibody against S100A12 (EN-RAGE, calgranulin-C), validated for IHC-P and ELISA on human samples. Ideal for studying inflammatory processes and alarmin signaling. Molecular weight 10 kDa.

[Learn More](#)

Parameter	Value
Product Name	S100A12 Rabbit pAb
Remarks/Alias	p6; CAGC; CGRP; MRP6; CAAF1; MRP-6; ENRAGE; S100A12
Species	Human
Gene ID (Human)	6283
Gene ID	6283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100A12 (NP_005612.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P80511
Protein Weight	10kDa
Shipping	Ice bag