

Anti-S100A12/Cgrp Rabbit Mab For Wb, If/Icc, Elisa - P80511

Item Number: CM0005569



Introduction

High-quality rabbit monoclonal antibody targeting human S100A12/CGRP, validated for WB, IF/ICC, ELISA. Suitable for inflammation and immune response research. Recombinant immunogen, 10 kDa.

[Learn More](#)

Parameter	Value
Product Name	S100A12/CGRP Rabbit mAb
Remarks/Alias	p6; CAGC; CGRP; MRP6; CAAF1; MRP-6; ENRAGE; S100A12/CGRP
Species	Human
Gene ID (Human)	6283
Gene ID	6283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100A12/CGRP (NP_005612.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P80511
Protein Weight	10kDa
Shipping	Ice bag