

Anti-S100A16 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Fq6

Item Number: CM0023350



Introduction

Rabbit polyclonal antibody against human S100A16, validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Suitable for studies of adipogenesis, calcium signaling, and cancer biomarker research.

[Learn More](#)

Parameter	Value
Product Name	S100A16 Rabbit pAb
Remarks/Alias	AAG13; S100F; DT1P1A7; S100A16
Species	Human
Gene ID (Human)	140576
Gene ID	140576
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96FQ6
Protein Weight	12kDa
Shipping	Ice bag