


Anti-S100A3 Rabbit Polyclonal Antibody For Wb And Elisa - P33764

Item Number: CM0025029



Introduction

S100A3 Rabbit Polyclonal Antibody validated for WB and ELISA, reacting with Human, Mouse, Rat. Targets S100A3, a calcium/zinc-binding protein involved in hair cuticle differentiation. Ideal for skin biology and hair research.

[Learn More](#)

Parameter	Value
Product Name	S100A3 Rabbit pAb
Remarks/Alias	S100E; S100A3
Species	Human
GeneID	6274
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human S100A3 (NP_002951.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P33764
Protein Weight	12kDa
Shipping	Ice bag