


## Anti-S100A6 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06703

Item Number: CM0008293



### Introduction

Rabbit monoclonal antibody targeting human S100A6 (calcyclin), validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat, enabling calcium signaling and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                       |
|------------------|-------------------------------------------------------------------------------------------------------------|
| Product Name     | S100A6 Rabbit mAb                                                                                           |
| Remarks/Alias    | 2A9; PRA; 5B10; CABP; CACY; S10A6; S100A6                                                                   |
| Species          | Human                                                                                                       |
| GeneID           | 6277                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human S100A6 (P06703) |
| Source           | Rabbit                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                       |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                    |
| Cross Reactivity | Human, Mouse, Rat                                                                                           |
| SWISS            | P06703                                                                                                      |
| Protein Weight   | 10kDa                                                                                                       |
| Shipping         | Ice bag                                                                                                     |