


Anti-S100A6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P06703

Item Number: CM0024253



Introduction

Rabbit polyclonal antibody against human S100A6 (Calcyclin) validated for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for calcium signaling and cancer research. Ideal for IVD raw material sourcing.

[Learn More](#)

Parameter	Value
Product Name	S100A6 Rabbit pAb
Remarks/Alias	2A9; PRA; 5B10; CABP; CACY; S10A6; S100A6
Species	Human
GeneID (human)	6277
GeneID	6277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human S100A6 (NP_055439.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06703
Protein Weight	10kDa
Shipping	Ice bag