


Anti-S100A6 Rabbit Polyclonal Antibody (Ko Validated) For Wb, If/Icc, Elisa - P06703

Item Number: CM0028378



Introduction

Rabbit polyclonal antibody against human S100A6, validated by knockout, suitable for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Involved in calcium signaling, cytoskeleton reorganization, and cell motility.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] S100A6 Rabbit pAb
Remarks/Alias	2A9; PRA; 5B10; CABP; CACY; S10A6; A6
Species	Human
GeneID (Human)	6277
GeneID	6277
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human S100A6 (NP_055439.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06703
Protein Weight	10kDa
Shipping	Ice bag