

Anti-S100A7/Psoriasin Polyclonal Antibody For Wb, Ihc-P, Elisa - P31151

Item Number: CM0022969



Introduction

Rabbit polyclonal antibody against human S100A7/Psoriasin, validated in WB, IHC-P, ELISA. Recombinant immunogen (aa 2-101), detects 11kDa protein. Ideal for studying psoriatic epidermis and SCC.

[Learn More](#)

Parameter	Value
Product Name	S100A7/Psoriasin Rabbit pAb
Remarks/Alias	PSOR1; S100A7c; S100A7/Psoriasin
Species	Human
Gene ID	6278 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-101 of human S100A7/Psoriasin (NP_002954.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P31151
Protein Weight	11kDa
Shipping	Ice bag