

Anti-S100A7/Psoriasin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P31151

Item Number: CM0006722



Introduction

Recombinant rabbit monoclonal antibody against human S100A7/Psoriasin (SWISS: P31151). Validated for WB, IHC-P, ELISA. Ideal for psoriasis and squamous cell carcinoma research.

[Learn More](#)

Parameter	Value
Product Name	S100A7/Psoriasin Rabbit mAb
Remarks/Alias	PSOR1; S100A7c; S100A7/Psoriasin
Species	Human
Gene ID (Human)	6278
Gene ID	6278
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P31151
Protein Weight	11 kDa
Shipping	Ice bag