

Anti-S100A7/Psoriasin Rabbit Polyclonal Antibody For IHC-P, Elisa - P31151

Item Number: CM0024150



Introduction

Rabbit polyclonal antibody against human S100A7/Psoriasin, a calcium-binding protein highly upregulated in psoriatic epidermis. Validated for IHC-P and ELISA; cross-reacts with human, mouse, rat. ~11 kDa.

[Learn More](#)

Parameter	Value
Product Name	S100A7/Psoriasin Rabbit pAb
Remarks/Alias	PSOR1; S100A7c; S100A7/Psoriasin
Species	Human
GeneID (Human)	6278
GeneID	6278
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-101 of human S100A7/Psoriasin (NP_002954.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31151
Protein Weight	11kDa
Shipping	Ice bag