


Anti-S100A8 Polyclonal Antibody For Wb, Ihc-P, Elisa - P05109

Item Number: CM0017546



Introduction

Rabbit polyclonal antibody targeting human S100A8, suitable for WB, IF/ICC, IHC-P, and ELISA. Shows cross-reactivity with mouse and rat. Ideal for inflammation, immunity, and cancer studies.

[Learn More](#)

Parameter	Value
Product Name	S100A8 Rabbit pAb
Remarks/Alias	P8; MIF; NIF; CAGA; CFAG; CGLA; L1Ag; MRP8; CP-10; MA387; 60B8AG; S100-A8; S100A8
Species	Human
GeneID (human)	6279
GeneID	6279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human S100A8 (NP_002955.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05109
Protein Weight	11kDa
Shipping	Ice bag