


Anti-S100A8 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05109

Item Number: CM0003026



Introduction

Rabbit monoclonal S100A8 antibody for WB, IHC-P, ELISA. Detects human, mouse, rat S100A8 (MRP8, Calgranulin-A), a key inflammatory mediator. Ideal for immunology, cancer research.

[Learn More](#)

Parameter	Value
Product Name	S100A8 Rabbit mAb
Remarks/Alias	P8; MIF; NIF; CAGA; CFAG; CGLA; L1Ag; MRP8; CP-10; MA387; 60B8AG; S100-A8; S100A8
Species	Human
Gene ID (Human)	6279
Gene ID	6279
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-93 of human S100A8 (P05109).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05109
Protein Weight	11 kDa
Shipping	Ice bag