


Anti-S100A8 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P05109

Item Number: CM0019461



Introduction

This S100A8 rabbit polyclonal antibody recognizes S100A8 (calgranulin-A/MRP8) in human, mouse, and rat. Suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Perfect for studying inflammatory processes, immune response, and cancer biomarkers.

[Learn More](#)

Parameter	Value
Product Name	S100A8 Rabbit pAb
Remarks/Alias	P8, MIF, NIF, CAGA, CFAG, CGLA, L1Ag, MRP8, CP-10, MA387, 60B8AG, S100-A8, S100A8
Species	Human
GeneID (Human)	6279
GeneID	6279
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human S100A8 (NP_002955.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05109
Protein Weight	11kDa
Shipping	Ice bag