

## Anti-S100A8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P27005

Item Number: CM0022011

### Introduction

Rabbit polyclonal antibody against mouse S100A8 (MRP8, calgranulin-A). Validated for Western blot, IHC-P, and ELISA. Detects mouse and rat S100A8, a 10 kDa inflammatory mediator.

[Learn More](#)

| Parameter        | Value                                             |
|------------------|---------------------------------------------------|
| Product Name     | S100A8 Rabbit pAb                                 |
| Remarks/Alias    | p8; B8Ag; CFAg; Caga; MRP8; CP-10; 60B8Ag; S100A8 |
| Species          | Mouse                                             |
| Gene ID (Human)  | 6279                                              |
| Gene ID          | 20201                                             |
| Immunogen        | Recombinant protein                               |
| Source           | Rabbit                                            |
| Category         | Polyclonal Antibodies                             |
| Application      | WB, IHC-P, ELISA                                  |
| Cross Reactivity | Mouse, Rat                                        |
| SWISS            | P27005                                            |
| Protein Weight   | 10kDa                                             |
| Shipping         | Ice bag                                           |