


Anti-S100A8 Rabbit Polymab For Wb, If-P, Ihc-P, Elisa - P05109 / P27005

Item Number: CM0009719



Introduction

S100A8 Rabbit PolymAb is a monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA in human and mouse. It targets S100A8, a calcium-binding protein involved in regulating inflammatory processes and immune response.

[Learn More](#)

Parameter	Value
Product Name	S100A8 Rabbit PolymAb®
Remarks/Alias	P8; MIF; NIF; CAGA; CFAG; CGLA; L1Ag; MRP8; CP-10; MA387; 60B8AG; S100-A8; S100A8; B8Ag; CFAg; Caga; 60B8Ag
Species	Human/Mouse
GeneID (Human)	6279
GeneID	6279/20201
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-93 of human/mouse S100A8 (P05109/P27005).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05109/P27005
Protein Weight	10 kDa / 11 kDa
Shipping	Ice bag