


Anti-S100A9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P31725

Item Number: CM0018762



Introduction

S100A9 rabbit polyclonal antibody detects calgranulin-B (MRP-14, p14) in mouse and human samples. Suitable for WB, IHC-P, ELISA. Target is a key regulator of inflammation and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	S100a9 Rabbit pAb
Remarks/Alias	p14; Cagb; GAGB; L1Ag; BEE22; MRP14; 60B8Ag
Species	Mouse
GeneID (Human)	6280
GeneID (Mouse)	20202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of mouse S100a9 (NP_033140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P31725
Protein Weight	13kDa
Shipping	Ice bag