

Anti-S100A9 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06702

Item Number: CM0030499



Introduction

High-quality rabbit polyclonal antibody against human S100A9 (calgranulin-B, MRP14), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for inflammation and cancer research.

[Learn More](#)

Parameter	Details
Product Name	S100A9 Rabbit pAb
Remarks/Alias	MIF; NIF; P14; CAGB; CFAG; CGLB; L1AG; LIAG; MRP14; 60B8AG; MAC387; S100-A9; S100A9
Species	Human
GeneID (human)	6280
GeneID	6280
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 20-114 of human S100A9 (NP_002956.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06702
Protein Weight	13kDa
Shipping	Ice bag