


Anti-S100B Rabbit Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P04271

Item Number: CM0018089



Introduction

S100B Rabbit pAb is a polyclonal antibody derived from rabbit, targeting human S100B protein. Recommended for WB, IP, IF-P, IHC-P, ELISA; cross-reacts with human, mouse, rat. Protein weight ~11 kDa.

[Learn More](#)

Parameter	Value
Product Name	S100B Rabbit pAb
Remarks/Alias	NEF; S100; S100-B; S100beta; S100B
Species	Human
GeneID (Human)	6285
GeneID	6285
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human S100B (NP_006263.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04271
Protein Weight	11kDa
Shipping	Ice bag