


Anti-S100P Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P25815

Item Number: CM0018422



Introduction

Rabbit polyclonal antibody against human S100P (MIG9), validated for Western blot, immunoprecipitation, immunofluorescence, immunohistochemistry, and ELISA. Excellent for studies on calcium signaling, microvilli biogenesis, and pancreatic ductal adenocarcinoma.

[Learn More](#)

Parameter	Value
Product Name	S100P Rabbit pAb
Remarks/Alias	MIG9; S100P
Species	Human
GeneID (Human)	6286
GeneID	6286
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P25815
Protein Weight	10kDa
Shipping	Ice bag