


Anti-S100P Rabbit Monoclonal Antibody For Wb, If/Icc, And Elisa - P25815

Item Number: CM0002054



Introduction

Rabbit monoclonal antibody targeting human S100P protein, validated for Western blot, immunofluorescence/immunocytochemistry, and ELISA. Recombinant immunogen; detects S100P in calcium signaling and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	S100P Rabbit mAb
Remarks/Alias	MIG9; S100P
Species	Human
GeneID (Human)	6286
GeneID	6286
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human S100P (NP_005971.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P25815
Protein Weight	10kDa
Shipping	Ice bag