

Anti-S100/S100A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23297

Item Number: CM0007012



Introduction

This S100/S100A1 Rabbit Monoclonal Antibody is validated for WB, IHC-P, IF/ICC, and ELISA. It detects human, mouse, and rat S100A1, a calcium-binding protein crucial for cardiomyocyte regulation and Ca²⁺ homeostasis.

[Learn More](#)

Parameter	Value
Product Name	S100/S100A1 Rabbit mAb
Remarks/Alias	S100; S100A; S100-alpha; S100/S100A1
Species	Human
GeneID (Human)	6271
GeneID	6271
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human S100/S100A1. (NP_006262.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23297
Protein Weight	11kDa
Shipping	Ice bag