

Anti-S100/S100A1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P23297

Item Number: CM0016294



Introduction

Rabbit polyclonal antibody against human S100/S100A1 for WB, IHC-P, ELISA. Recognizes human, mouse, rat. S100A1 is a cardiac calcium-binding protein regulating Ca^{2+} homeostasis.

[Learn More](#)

Parameter	Value
Product Name	S100/S100A1 Rabbit pAb
Remarks/Alias	S100; S100A; S100-alpha; S100/S100A1
Species	Human
Gene ID (Human)	6271
Gene ID	6271
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-94 of human S100/S100A1 (NP_006262.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23297
Protein Weight	11kDa
Shipping	Ice bag