

Anti-S1Pr2 Rabbit Polyclonal Antibody For Wb And Elisa - O95136

Item Number: CM0019536



Introduction

High-quality rabbit polyclonal antibody targeting human S1PR2 (O95136), validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for S1P signaling, GPCR, and hearing loss research.

[Learn More](#)

Parameter	Value
Product Name	S1PR2 Rabbit pAb
Remarks/Alias	EDG5; H218; LPB2; S1P2; AGR16; EDG-5; DFNB68; Gpcr13; S1PR2
Species	Human
GeneID (Human)	9294
GeneID	9294
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 250-353 of human S1PR2 (NP_004221.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95136
Protein Weight	39kDa
Shipping	Ice bag