


Anti-S1Pr3 Rabbit Monoclonal Antibody For Wb, Elisa - Q99500

Item Number: CM0002339



Introduction

High-quality anti-S1PR3 rabbit monoclonal antibody validated for WB and ELISA applications. Reacts with human, mouse, and rat S1PR3 (EDG3). Ideal for studying sphingosine-1-phosphate receptor signaling and related biological processes.

[Learn More](#)

Parameter	Value
Product Name	S1PR3 Rabbit mAb
Remarks/Alias	EDG3; LPB3; S1P3; EDG-3; C9orf47; C9orf108; bA791021.3; S1PR3
Species	Human
GeneID (Human)	1903
GeneID	1903
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human S1PR3 (Q99500).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99500
Protein Weight	42kDa
Shipping	Ice bag