


Anti-S1Pr5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H228

Item Number: CM0018987



Introduction

S1PR5 rabbit polyclonal antibody for WB and ELISA. Targets human sphingosine 1-phosphate receptor 5 (EDG8), a GPCR involved in diverse physiological processes. For neuroscience and immunology studies.

[Learn More](#)

Parameter	Value
Product Name	S1PR5 Rabbit pAb
Remarks/Alias	EDG8; S1P5; Edg-8; SPPR-1; SPPR-2; S1PR5
Species	Human
Gene ID	53637
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 307-398 of human S1PR5 (NP_110387.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H228
Protein Weight	42kDa
Shipping	Ice bag