

Anti-S6 Ribosomal Protein (Rps6) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62753

Item Number: CM0002765



Introduction

Rabbit monoclonal anti-RPS6 antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets 29 kDa 40S ribosomal protein eS6. Ideal for translation and cell growth studies.

[Learn More](#)

Parameter	Value
Product Name	S6 Ribosomal Protein (RPS6) Rabbit mAb
Remarks/Alias	S6; eS6; S6 Ribosomal Protein (RPS6)
Species	Human
Gene ID (Human)	6194
Gene ID	6194
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RPS6 (NP_001001.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag