


Anti-S6 Ribosomal Protein (Rps6) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62753

Item Number: CM0029367



Introduction

Rabbit polyclonal antibody against human S6 Ribosomal Protein (RPS6). Suitable for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. For ribosome and cell growth studies.

[Learn More](#)

Parameter	Value
Product Name	S6 Ribosomal Protein (RPS6) Rabbit pAb
Remarks / Alias	S6; eS6; S6 Ribosomal Protein (RPS6)
Species	Human
Gene ID (Human)	6194
Gene ID	6194
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-130 of human RPS6 (NP_001001.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62753
Protein Weight	29kDa
Shipping	Ice bag