


Anti-Saa1 Polyclonal Antibody For Wb And Elisa - P0Dji8

Item Number: CM0015370



Introduction

Rabbit polyclonal antibody targeting human SAA1 (Serum Amyloid A1). Validated for WB and ELISA applications. Cross-reacts with mouse and rat. Ideal for acute phase protein and amyloidosis research.

[Learn More](#)

Parameter	Value
Product Name	SAA1 Rabbit pAb
Remarks/Alias	SAA; PIG4; SAA2; TP53I4; SAA1
Species	Human
GeneID (human)	6288
GeneID	6288
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SAA1 (NP_000322.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P0DJI8
Protein Weight	12 kDa
Shipping	Ice bag