

Anti-Saal1 Polyclonal Antibody For Wb, Elisa - Q96Er3

Item Number: CM0014377

Introduction

Rabbit polyclonal antibody targeting human SAAL1 (SPACIA1), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for research on synovial fibroblast proliferation and inflammatory arthritis.

[Learn More](#)

Parameter	Value
Product Name	SAAL1 Rabbit pAb
Remarks/Alias	SPACIA1; SAAL1
Species	Human
GeneID (Human)	113174
GeneID	113174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-474 of human SAAL1 (NP_612430.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96ER3
Protein Weight	54kDa
Shipping	Ice bag