


Anti-Sae1 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Ube0

Item Number: CM0003310



Introduction

SAE1 Rabbit mAb validated for WB, ELISA. Targets human SAE1 (SUMO-activating enzyme subunit 1), a key E1 enzyme in SUMOylation, essential for nuclear protein regulation and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	SAE1 Rabbit mAb
Remarks/Alias	AOS1, SUA1, UBLE1A, HSPC140, SAE1
Species	Human
GeneID (Human)	10055
GeneID	10055
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SAE1 (Q9UBE0).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UBE0
Protein Weight	38kDa
Shipping	Ice bag