


Anti-Sae1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ube0

Item Number: CM0030252



Introduction

Rabbit polyclonal antibody targeting human SAE1 (SUMO-activating enzyme subunit 1). Suitable for Western blot and ELISA, detects human and mouse SAE1. For SUMOylation and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	SAE1 Rabbit pAb
Remarks/Alias	AOS1; SUA1; UBLE1A; HSPC140; SAE1
Species	Human
GeneID (Human)	10055
GeneID	10055
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-346 of human SAE1 (NP_005491.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UBE0
Protein Weight	38kDa
Shipping	Ice bag