

Anti-Safb Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15424

Item Number: CM0007022

Introduction

Anti-SAFB rabbit monoclonal antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Recombinant immunogen for transcription regulation and estrogen receptor corepressor research.

[Learn More](#)

Parameter	Value
Product Name	SAFB Rabbit mAb
Remarks/Alias	HAP, HET, SAF-B, SAFB1, SAB-B1, SAF-B1, SAFB
Species	Human
GeneID (Human)	6294
GeneID	6294
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 815-915 of human SAFB (NP_002958.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15424
Protein Weight	103kDa
Shipping	Ice bag