

## Anti-Safb2 Rabbit Polyclonal Antibody For Wb, Elisa - Q14151

Item Number: CM0027038

### Introduction

Anti-SAFB2 rabbit polyclonal antibody validated for Western blot and ELISA. Detects SAFB2 (~107 kDa) in human, mouse, and rat samples. Targets a scaffold attachment factor B2 involved in chromatin organization and transcription regulation. Ideal for cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | SAFB2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | SAFB2                                                                                                            |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 9667                                                                                                             |
| GeneID           | 9667                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human SAFB2 (NP_055464.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | Q14151                                                                                                           |
| Protein Weight   | 107kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                          |