

Anti-Safb2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14151

Item Number: CM0015301

Introduction

SAFB2 Rabbit pAb, a polyclonal antibody for WB, IHC-P, ELISA, detects human, mouse, and rat SAFB2, a scaffold attachment factor involved in estrogen receptor repression and cell proliferation inhibition.

[Learn More](#)

Parameter	Value
Product Name	SAFB2 Rabbit pAb
Remarks/Alias	SAFB2
Species	Human
GeneID (Human)	9667
GeneID	9667
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human SAFB2 (NP_055464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14151
Protein Weight	107kDa
Shipping	Ice bag