

Anti-Sag (S-Arrestin) Polyclonal Antibody - P10523

Item Number: CM0016965

Introduction

Rabbit polyclonal antibody against human SAG (S-arrestin) for Western blot (WB), immunohistochemistry (IHC-P), and ELISA. Cross-reacts with mouse and rat. Suitable for rhodopsin signaling and retinal research.

[Learn More](#)

Parameter	Value
Product Name	SAG Rabbit pAb
Remarks/Alias	RP47; RP96; S-AG; SAG
Species	Human
GeneID	6295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-405 of human SAG (NP_000532.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10523
Protein Weight	45kDa
Shipping	Ice bag