

Anti-Sall4 Monoclonal Antibody For Ihc-P And Elisa - Q9Ujq4

Item Number: CM0005886



Introduction

Rabbit monoclonal antibody targeting human Sall4 for IHC-P and ELISA applications. Ideal for stem cell maintenance and leukemia studies. Cross-reacts with human. ~112 kDa.

[Learn More](#)

Parameter	Value
Product Name	Sall4 Rabbit mAb
Remarks/Alias	DRRS; IVIC; HSAL4; ZNF797; Sall4
Species	Human
GeneID (human)	57167
Immunogen	Recombinant protein of human Sall4
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9UJQ4
Protein Weight	112kDa
Shipping	Ice bag