

Anti-Sall4 Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Ujq4

Item Number: CM0018839



Introduction

This anti-Sall4 rabbit polyclonal antibody targets human Sall4 (UniProt Q9UJQ4), validated for WB, IF-P, IHC-P, and ELISA, with cross-reactivity to mouse and rat, supporting stem cell research.

[Learn More](#)

Parameter	Value
Product Name	Sall4 Rabbit pAb
Remarks/Alias	DRRS; IVIC; HSAL4; ZNF797; Sall4
Species	Human
GeneID (Human)	57167
GeneID	57167
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 853-1053 of human Sall4 (NP_065169.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UJQ4
Protein Weight	112kDa
Shipping	Ice bag