

Anti-Samhd1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y3Z3

Item Number: CM0006074



Introduction

SAMHD1 Rabbit mAb validated for WB, IHC-P, ELISA with human cross reactivity. Targets deoxynucleoside triphosphate triphosphohydrolase involved in antiviral defense and DNA repair; suitable for immunology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	SAMHD1 Rabbit mAb
Remarks/Alias	DCIP; CHBL2; HDDC1; MOP-5; SBB188; hSAMHD1; SAMHD1
Species	Human
Gene ID	25939
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 527-626 of human SAMHD1 (NP_056289.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9Y3Z3
Protein Weight	72kDa
Shipping	Ice bag